



RESEARCH KITS

ELISA AND RIA

For Small Sample Volumes

ELISA TESTS INCLUDE:

- ADRENALINE
- CITRULLINE
- DOPAMINE
- GABA
- GLUTAMATE
- GLYCINE
- HISTAMINE
- KYNURENINE
- NORADRENALINE
- SEROTONIN
- TRYPTOPHAN

RIA TEST INCLUDE:

- MELATONIN

FEATURES

- High sensitivity
- Small sample volumes
- Flexible protocols
- Designed to test various sample types from multiple species
- Customer Service! – We use the kits we sell, and we always offer a 24 hour turn-around on almost all technical questions



IMMUNOASSAYS AND SERVICES

BIOGENIC AMINES & NEUROSCIENCE | ENDOCRINOLOGY | FOOD SAFETY

RESEARCH KITS

Elisa & RIA

CAT NO.	PRODUCT NAME	SAMPLE TYPE AND VOLUME	SENSITIVITY
BA E-5100	Adrenaline Research ELISA	1 - 750 µl, all biological fluids	3.3 pg/ml (depending on sample volume)
BA E-5200	Noradrenaline Research ELISA	1 - 750 µl, all biological fluids	1.3 pg/ml (depending on sample volume)
BA E-5300	Dopamine Research ELISA	1 - 750 µl, all biological fluids	3.3 pg/ml (depending on sample volume)
BA E-5400	2-CAT (A-N) Research ELISA	1 - 750 µl, all biological fluids	Adrenaline: 6.6 pg/ml Noradrenaline: 2.6 pg/ml (depending on sample volume)
BA E-5500	2-CAT (N-D) Research ELISA	1 - 750 µl, all biological fluids	Noradrenaline 2.6 pg/ml Dopamine 6.6 pg/ml (depending on sample volume)
BA E-5600	3-CAT Research ELISA	1 - 750 µl, all biological fluids	Adrenaline: 10 pg/ml plasma Noradrenaline: 4 pg/ml plasma Dopamine: 10 pg/ml plasma (depending on sample volume)
BA E-2800	Citrulline ELISA	50 µl serum, plasma, urine other biological fluids	0.23 µg/ml
BA E-2500	GABA ELISA	100 µl urine, 300 µl serum, plasma, cell cultures and tissue homogenates	49 ng/ml (urine samples)
BA E-2400	Glutamate ELISA	100 µl urine, 100 µl serum, plasma, other biological fluids	0.3 µg/ml
BA E-2100	Glycine ELISA	20 µl urine, other biological fluids	3.3 µg/ml
BA E-1700	Histamine ELISA ^{Fast Track}	20 µl cell culture, whole blood, or urine	0.2 ng/ml
BA E-5800	Histamine Research ELISA	10 µl, all biological fluids	0.2 ng/ml
BA E-2200	Kynurenine ELISA	20 µl serum or plasma, other biological fluids	45.7 ng/ml
BA R-3900	Melatonin Research RIA	25-200 µl, all biological fluids	5.9 pg/ml (depending on sample volume)
BA E-5900	Serotonin Research ELISA	1-100 µl, all biological fluids	6,2 pg/ml (depending on sample volume)
BA E-2700	Tryptophan ELISA	20 µl urine, serum, plasma, other biological fluids	1.2 µg/ml