

New Arrival

Elabscience®

Affinity Agarose & Immunomagnetic Beads

Best Tools for IP/CoIP!

10-Years of Research and Development Fully Verified Qualified and Reliable



3 Easy Steps to Finish Your Immunoprecipitation




Fig1. Affinity Agarose Operation Process Chart




Fig2. Immunomagnetic Beads Operation Process Chart

Why choose our IP/CoIP products?

★High Affinity ★High Sensitivity ★High Specificity




Fig3. GFP-DYKDDDDK Protein Purified by DYKDDDDK Affinity Agarose

A crude cell extract was run over the DYKDDDDK affinity agarose, only the protein of interest was bound, and then eluted by 3x DYKDDDDK peptide.

Over 10 Categories

Most commonly-used tags are covered.

Best Services

Technique support full of experience and efficiency.

Original Manufacturer

Sufficient supplies.

» Affinity Agarose & Immunomagnetic Beads IP/CoIP Kits

- For IP/CoIP assays.
- Containing all of the reagents needed, to assure your result.

Cat. No.	Product Name	Size/Price
EA-IP-K001	Anti-DYKDDDDK-Tag FAST IP Kit	25 Test
EA-IP-K002	Anti-HA Tag (YPYDVPDYA) FAST IP Kit	25 Test
EA-IP-K003	Anti-MYC Tag (EQKLISEEDL) FAST IP Kit	25 Test
EA-IP-K008	GST Pull-Down Kit(Agarose)	25 Test
EA-IP-K007	Protein A/G Affinity Agarose	50 Test
EA-IP-K007M	Protein A/G Magnetic Bead	50 Test

» Affinity Agarose & Immunomagnetic Beads-Reagents

- For IP/CoIP assays, and tagged protein purification as well.
- Made by covalent coupling of high-quality Tag antibody and Sepharose 4B agarose granules., with the characteristics of high binding capacity of protein, high specificity, stability and repeatability.

Cat. No.	Product Name	Size/Price
EA-IP-001	Anti-DYKDDDDK Tag (DYKDDDDK) Affinity Agarose	2mL
EA-IP-002	Anti-HA Tag (YPYDVPDYA) Affinity Agarose	2mL
EA-IP-003	Anti-MYC Tag (EQKLISEEDL) Affinity Agarose	2mL
EA-IP-004	Anti-GFP Tag Affinity Agarose	2mL
EA-IP-004N	Anti-GFP Tag Nanobody Affinity Agarose	2mL
EA-IP-005N	Anti-mCherry Nanobody Affinity Agarose	2mL
EA-IP-006N	Anti-MBP Nanobody Affinity Agarose	2mL
EA-IP-001M	Anti-DYKDDDDK Tag (DYKDDDDK) Immunomagnetic Beads	1mL
EA-IP-002M	Anti-HA Tag (YPYDVPDYA) Immunomagnetic Beads	1mL
EA-IP-003M	Anti-MYC (EQKLISEEDL) Immunomagnetic Beads	1mL
EA-IP-004MN	Anti-GFP Nanobody Immunomagnetic Beads	1mL
EA-IP-005MN	Anti-mCherry Nanobody Immunomagnetic Beads	1mL
EA-IP-007	Protein A/G Agarose	2mL
EA-IP-007M	Protein A/G Magnetic Beads	1mL

» Tagged Protein Purification-Kits

- For tagged protein purification
- Containing all of the reagents needed, including affinity columns, binding, wash, elution and storage buffers.

Cat. No.	Chinese Name	Size/Price
EA-TP-K001	DYKDDDDK-tagged Protein Purification Kit	1 kit
EA-TP-K002	HA (YPYDVPDYA)-tagged Protein Purification Kit	1 kit
EA-TP-K003	MYC (EQKLISEEDL)-tagged Protein Purification Kit	1 kit

Elabscience Biotechnology Inc.

Tel: 1-240-252-7368
Fax: 1-240-252-7376

Web: www.elabscience.com
Email: techsupport@elabscience.com